

28-N1 34 154 34 23 905

Delivery System Usable Length & Device Category

Delivery System Outer Diameter Fr

Stent Distal Diameter

Stent Covered Length

Stent Proximal Diameter

Device Type (N=Non Bare Stent) & Modification No.

Product Identification

DS = Delivery System • Fr = French Size • O.D. = Outer Diameter

RELAY™ NBS
THORACIC STENT-GRAFT

Non Bare Stent

- New Non Bare Stent configuration with **covered proximal end**
- Wide treatment range
- Controlled and precise stent graft deployment

When treating Aneurysms and Dissections in Thoracic Aorta
RELY on RELAY

Werfen Group | Bolton Medical


This product is certified by the TÜV SÜD product SERVICE (0123) according to 93/42 EEC.
Manufacturer: Bolton Medical España S.L.U. C/ Newton, 18-24,
08635 Sant Esteve de Sesrovières, Barcelona, Spain. +34- 93-817-63-00

Relay NBS not for sale in the U.S.



Bolton Medical

Dimension and Reference Guide



RELAY™ NBS
THORACIC STENT-GRAFT

RELY on RELAY

28-MK0021_04_09

Dimension and Reference Guide

STRAIGHT GRAFT		STENT GRAFT			DELIVERY SYSTEM			VESSEL
Catalog Number		Proximal/Distal Diameter (mm)	Length Category (mm)	Covered Length (mm)	DS.O.D. (Fr)	Usable Length (cm)	Tip to Graft (cm)	Vessel Diameter (mm)
28-N1 26 104 26 22 90S	26	100	104	22	22	90	4.9	22-23
28-N1 26 164 26 22 90S	26	150	164	22	22	90	4.9	22-23
28-N1 26 204 26 22 90S	26	200	204	22	22	90	4.9	22-23
28-N1 26 259 26 22 90S	26	250	259	22	22	90	4.9	22-23
28-N1 28 104 28 22 90S	28	100	104	22	22	90	4.9	24-25
28-N1 28 164 28 22 90S	28	150	164	22	22	90	4.9	24-25
28-N1 28 204 28 22 90S	28	200	204	22	22	90	4.9	24-25
28-N1 28 259 28 22 90S	28	250	259	22	23	90	5.1	24-25
28-N1 30 104 30 22 90S	30	100	104	22	22	90	4.9	26-27
28-N1 30 164 30 22 90S	30	150	164	22	22	90	4.9	26-27
28-N1 30 209 30 23 90S	30	200	209	23	23	90	5.1	26-27
28-N1 30 259 30 23 90S	30	250	259	23	23	90	5.1	26-27
28-N1 32 104 32 23 90S	32	100	104	23	23	90	5.1	28-29
28-N1 32 164 32 23 90S	32	150	164	23	23	90	5.1	28-29
28-N1 32 209 32 23 90S	32	200	209	23	23	90	5.1	28-29
28-N1 32 259 32 24 90S	32	250	259	24	24	90	5.3	28-29
28-N1 34 109 34 23 90S	34	100	109	23	23	90	5.1	30-31
28-N1 34 154 34 23 90S	34	150	154	23	23	90	5.1	30-31
28-N1 34 209 34 24 90S	34	200	209	24	24	90	5.3	30-31
28-N1 34 259 34 24 90S	34	250	259	24	24	90	5.3	30-31
28-N1 36 109 36 24 90S	36	100	109	24	24	90	5.3	32-33
28-N1 36 154 36 24 90S	36	150	154	24	24	90	5.3	32-33
28-N1 36 199 36 24 90S	36	200	199	24	24	90	5.3	32-33
28-N1 36 259 36 25 90S	36	250	259	25	25	90	5.5	32-33
28-N1 38 109 38 24 90S	38	100	109	24	24	90	5.3	34
28-N1 38 154 38 24 90S	38	150	154	24	24	90	5.3	34
28-N1 38 199 38 25 90S	38	200	199	25	25	90	5.5	34
28-N1 38 259 38 25 90S	38	250	259	25	25	90	5.5	34
28-N1 40 114 40 25 90S	40	100	114	25	25	90	5.5	35-36
28-N1 40 154 40 25 90S	40	150	154	25	25	90	5.5	35-36
28-N1 40 204 40 25 90S	40	200	204	25	25	90	5.5	35-36
28-N1 40 259 40 25 90S	40	250	259	25	25	90	5.5	35-36
28-N1 42 114 42 25 90S	42	100	114	25	25	90	5.5	37-38
28-N1 42 159 42 25 90S	42	150	159	25	25	90	5.5	37-38
28-N1 42 204 42 25 90S	42	200	204	25	25	90	5.5	37-38
28-N1 42 259 42 25 90S	42	250	259	25	25	90	5.5	37-38
28-N1 44 114 44 25 90S	44	100	114	25	25	90	5.5	39-40
28-N1 44 164 44 25 90S	44	150	164	25	25	90	5.5	39-40
28-N1 44 209 44 25 90S	44	200	209	25	25	90	5.5	39-40
28-N1 44 259 44 25 90S	44	250	259	25	25	90	5.5	39-40
28-N1 46 114 46 26 90S	46	100	114	26	26	90	5.5	41-42
28-N1 46 164 46 26 90S	46	150	164	26	26	90	5.5	41-42
28-N1 46 209 46 26 90S	46	200	209	26	26	90	5.5	41-42
28-N1 46 259 46 26 90S	46	250	259	26	26	90	5.5	41-42

Note: Length 100 and 250 available on demand (delivery time up to 3 weeks)

TAPERED/CONICAL GRAFT		STENT GRAFT			DELIVERY SYSTEM			VESSEL
Catalog Number		Proximal/Distal Diameter (mm)	Length Category (mm)	Covered Length (mm)	DS.O.D. (Fr)	Usable Length (cm)	Tip to Graft (cm)	Vessel Diameter (mm)
28-N1 34 154 30 23 90S	34/30	150	154	23	23	90	5.1	30-31
28-N1 34 209 30 24 90S	34/30	200	209	24	24	90	5.3	30-31
28-N1 36 154 32 24 90S	36/32	150	154	24	24	90	5.3	32-33
28-N1 36 199 32 24 90S	36/32	200	199	24	24	90	5.3	32-33
28-N1 38 154 34 24 90S	38/34	150	154	24	24	90	5.3	34
28-N1 38 199 34 25 90S	38/34	200	199	25	25	90	5.5	34
28-N1 40 154 36 25 90S	40/36	150	154	25	25	90	5.5	35-36
28-N1 40 204 36 25 90S	40/36	200	204	25	25	90	5.5	35-36
28-N1 42 159 38 25 90S	42/38	150	159	25	25	90	5.5	37-38
28-N1 42 204 38 25 90S	42/38	200	204	25	25	90	5.5	37-38
28-N1 44 164 40 25 90S	44/40	150	164	25	25	90	5.5	39-40
28-N1 44 209 40 25 90S	44/40	200	209	25	25	90	5.5	39-40
28-N1 46 164 42 26 90S	46/42	150	164	26	26	90	5.5	41-42
28-N1 46 209 42 26 90S	46/42	200	209	26	26	90	5.5	41-42